



CERTIFICATE OF MATERIAL ORIGIN

Date: 10/27/2017

Revision: 02

Product: Poly-L-Lysine, Solution, 0.1 mg/ml

Supplier	Advanced BioMatrix, Inc.
Supplier Catalog No.	5048
Supplier Location	California, USA
Name/Grade	Poly-L-Lysine, Solution, 0.1 mg/ml
Animal Origin (including human)	
-Animal Species	Not applicable
-Tissue Origin	Not applicable
-County of Origin	Not applicable
Non-animal Source of Origin	
-Country of Origin	United States
Animal Origin Material Adventitious Agent Risk Reduction Steps	Not applicable
Method of Manufacturing	Synthetic Processing
Briefly describe Method of Manufacturing	Poly-L-lysine is a positively charged amino acid polymer with approximately one HBr per lysine residue. The hydrobromide allows the poly-L-lysine to be in a crystalline form soluble in water. A small amount of product may be found in the beta structure because the HBr interferes with hydrogen bonding between amino and either the carboxyl groups or N or O containing moieties.

The above information has been completed by an authorized representative of management and is complete and correct to my knowledge.

President

October 27, 2017

Corporate Headquarters
P.O. Box 502403
San Diego, CA 92150-2403

(800) 883-8220 (760) 929-0755
(510) 217-3452 Fax
www.advancedbiomatrix.com

Manufacturing Facility
5930 Sea Lion Place, Ste. 100
Carlsbad, CA 92010